


MACHEREY-NAGEL

NANOCOLOR[®] Robot tests

Water Analysis



NEW!
COD up to
15 g/L O₂

Skalar

NANOCOLOR[®] – Automated water analysis

- NANOCOLOR[®] tube tests on the SKALAR Robot SP2000^{series}
- Fully automatic photometric water analysis including sample preparation and digestion
- Analysis of all important waste water parameters

MACHEREY-NAGEL

www.mn-net.com



NANOCOLOR® Robot tests

Automated water analysis with NANOCOLOR®

The companies MACHEREY-NAGEL and Skalar joined forces and designed the SP2000^{series}, a robotic analyzer for the fully automated water analysis in the laboratory. This water analyzer processes a variety of NANOCOLOR® tests and automates all the necessary handling steps of the photometric test kit procedure such as sample pipetting, (de)-capping, addition of reagents, mixing, heating, cooling and measurement. The flexibility and versatility of the test kit platform is for example reflected in the possibility to process from 24 to 336 test tubes at once by using different test tube racks. The instrument can be configured for the analysis of one test kit application or for multiple test kit applications per sequence such as combinations of COD, total phosphate, total nitrogen and other parameters.



Fully automatic analyzable NANOCOLOR® tests

Test	REF	Measuring range	No. of tests
Ammonium 3	985603	0.04–2.30 mg/L NH ₄ -N	20
Ammonium 50	985605	1–40 mg/L NH ₄ -N	20
Chloride 50 *	985621	0.5–50.0 mg/L Cl ⁻	20
Chloride 200 *	985619	5–200 mg/L Cl ⁻	20
COD 60 *	ISO 15705 985622	5–60 mg/L O ₂	20
COD 160 *	ISO 15705 985626	10–160 mg/L O ₂	20
COD 600 *	ISO 15705 985630	50–600 mg/L O ₂	20
COD 1500 *	ISO 15705 985629	100–1500 mg/L O ₂	20
COD 15000 *	985628	1.0–15.0 g/L O ₂	20
COD LR 150 *	ISO 15705 985636	3–150 mg/L O ₂	19
COD HR 1500 *	ISO 15705 985638	20–1500 mg/L O ₂	19
Nitrate 8 *	985665	0.30–8.00 mg/L NO ₃ -N	20
Nitrate 50 *	985664	0.3–22.0 mg/L NO ₃ -N	20
Nitrite 2 *	985668	0.003–0.460 mg/L NO ₂ -N	20
Nitrite 4	985669	0.1–4.0 mg/L NO ₂ -N	20
Phenolic index 5 *	985674	0.2–5.0 mg/L Phenol	20
total Phosphate 1 *	985676	0.05–1.50 mg/L P	20
total Phosphate 5 *	985681	0.20–5.00 mg/L P	20
total Phosphate 15 *	985680	0.30–15.00 mg/L P	20
ortho Phosphate 1 *	985607	0.05–1.50 mg/L P	20
ortho Phosphate 15 *	985657	0.30–15.00 mg/L P	20
total Nitrogen TN ₀ 22 *	985683	0.5–22.0 mg/L N	20
total Nitrogen TN ₀ 220 *	985688	5–220 mg/L N	20
Sulfide 3 *	985673	0.05–3.00 mg/L S ²⁻	20

Only with MN: COD up to 15 g/L O₂

* This product contains harmful substances which must be specially labelled as hazardous. For detailed information please see SDS.

Contact:

Technical Support Water Analysis:

+49 24 21 969-332 · csc@mn-net.com



www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

