

# NANOCOLOR® Ammonium 3

Amonio


690 nm

Method(e) / Método


0031	0.04 - 2.30 mg/l NH <sub>4</sub> -N
0032	0.05 - 3.00 mg/l NH <sub>4</sub> <sup>+</sup>
0033	0.05 - 3.00 mg/l NH <sub>3</sub>

**Test 0-03**


REF 985 003




**1**  
4.0 ml Probe  
Sample  
Enchantillon  
Muestra




**2**  
1 x NANOFIX R2




**3**  
Schütteln  
Shake  
Agiter  
Agitar




**4**  
Säubern  
Clean  
Nettoyer  
Limpiar



**5**  
15'00 min



**6**  
Messung  
Measurement  
Mesure  
Medición



1:1

Meerwasser  
Sea water  
Eau de mer  
Agua de mar