

NANOCOLOR® Ammonium 50

Amonio


690 nm

Method(e) / Método


0051	1 - 40 mg/l $\text{NH}_4\text{-N}$
0052	1 - 50 mg/l NH_4^+
0053	1 - 50 mg/l NH_3

Test 0-05


REF 985 005




200 µl Probe
Sample
Enchantillon
Muestra




1 x NANOFIX R2




Schütteln
Shake
Agiter
Agitar




Säubern
Clean
Nettoyer
Limpiar



15'00 min



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar