

NANOCOLOR® Kationische Tenside 4

Cationic surfactants / Tensio-actifs
cationiques / Tensioactivos catiónicos

620 nm


Method(e) / Método

0341


0.20 - 4.00 mg/l CTBA

Test 0-34


REF 985 034




1
4.0 ml Probe
Sample
Enchantillon
Muestra




2
500 µl R2




3
Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo




4
Säubern
Clean
Nettoyer
Limpiar



5
10'00 min
Phasentrennung abwarten
wait for phase separation
laisser séparer les phases
esperar la separación de fases



6
Messung
Measurement
Mesure
Medición



1:20
Meerwasser
Sea water
Eau de mer
Agua de mar