

NANOCOLOR® Formaldehyd(e) 8

Formaldehído

585 nm


Method(e) / Método

0411


0.1 - 8.0 mg/l HCHO

Test 0-41


REF 985 041




2.0 ml Probe
Sample
Enchantillon
Muestra




1 x R2




Sicherheitsgefäß
Safety bottle
Réipient de sécurité
Recipiente de seguridad



5'00 min



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar