

NANOCOLOR® Molybdän 40

Molybdenum / Molybdène / Molibdeno

345 nm


Method(e) / Método

0561 1.0 - 30.0 mg/l Mo(VI)


0562 1.6 - 50.0 mg/l MoO₄²⁻

Test 0-56


REF 985 056




1
4.0 ml Probe
Sample
Enchantillon
Muestra




2
Schütteln
Shake
Agiter
Agitar



3
Säubern
Clean
Nettoyer
Limpiar



4
5'00 min



5
Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar