

NANOCOLOR® Sauerstoff 12

Oxygen / Oxygène / Oxígeno

436 nm


Method(e) / Método

0821

0.5 - 12.0 mg/l O₂

Test 0-82


REF 985 082




Probe
Sample
Enchantillon
Muestra




Nullwert
Blank value
Zéro
Cero




2 x ◊ R1




2 x ◊ R2




Schütteln
Shake
Agiter
Agitar




2'00 min




5 x ◊ R3



Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición

Meerwasser
Sea water
Eau de mer
Agua de mar

