

NANOCOLOR® Thiocyanat 50

Thiocyanate / Tiocianato

470 nm


Method(e) / Método

0911


0.5 - 50.0 mg/l SCN⁻

Test 0-91


REF 985 091




4.0 ml Probe
Sample
Enchantillon
Muestra




Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición



1:1
Meerwasser
Sea water
Eau de mer
Agua de mar