

visicolor® ECO Sulfat

Sulphate / Sulfate / Sulfato

Test 5-92


REF 931 292

436 nm


Method(e) / Método

5921


20–200 mg/L SO_4^{2-}




Probe / sample
échantillon / muestra




Küvetten spülen
Rinse test tubes
Rincer les éprouvettes
Lavar los tubos




10 mL Probe
Sample
Echantillon
Muestra




10 $\text{SO}_4\text{-1}$




Schütteln
Shake
Agiter
Agitar




1 $\text{SO}_4\text{-2}$




Schütteln
Shake
Agiter
Agitar




1'00 min



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar