

# pH 6.5 – 8.2

# Test 0-72

REF 918 72



Rundküvettest / Tube test

Test en cuves rondes / Prueba en cubeta cilíndrica

436 / 540 nm


Method(e) / Método

0721


pH 6.5 - 8.2

Nullwert / Blanc value / Zéro / Cero


Messwert / Sample / Echantillon / Muestra




1  
10 ml Probe  
Sample  
Echantillon  
Muestra




1  
10 ml Probe  
Sample  
Echantillon  
Muestra




2  
1 x NANOFIX pH




3  
Säubern  
Clean  
Nettoyer  
Limpiar




4  
Schütteln  
Shake  
Agiter  
Agitar




5  
Säubern  
Clean  
Nettoyer  
Limpiar



6  
Messung  
Measurement  
Mesure  
Medición



7  
3'00 min



8  
Messung  
Measurement  
Mesure  
Medición