

CSB 160

COD/DCO/DQO

Test 0-26

REF 985 026



Rundküvettentest/ Tube test


Test en cuves rondes / Prueba en cubeta cilíndrica

15–160 mg/L


436 nm

Method(e) / Método


0261




2.0 mL Probe
Sample
Echantillon
Muestra




Sicherheitsgefäß
Safety bottle
Réceptient de sécurité
Recipiente de seguridad




148 °C / 2 h




Schütteln
Shake
Agiter
Agitar



Abkühlen lassen
Allow to cool
Laisser refroidir
Dejar enfriar



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición