

# organische Säuren 3000

## organic Acids / Acides organiques

### Acidos orgánicos

Rundküvettentest / Tube test

Test en cuves rondes / Prueba en cubeta cilíndrica

30–3000 mg/L  $\text{CH}_3\text{COOH}$   
0.5–50.0 mmol/L  $\text{CH}_3\text{COOH}$

**Test 0-50**

REF 985 050




470 nm


Method(e) / Método

0501


0502




500  $\mu\text{L}$  R2




1.0 mL Probe  
Sample  
Echantillon  
Muestra




Schütteln  
Shake  
Agiter  
Agitar




100 °C / 10 min




Sofort abkühlen  
Cool **immediately**  
Refródir **immédiatement**  
Enfriar **inmediata**




1 x R3



Schütteln  
Shake  
Agiter  
Agitar



1.0 mL R4



Schütteln  
Shake  
Agiter  
Agitar