

Molybdän 40

Molybdenum / Molybdène / Molibdeno

Test 0-56

REF 985 056




Rundküvettest / Tube test


Test en cuves rondes / Prueba en cubeta cilíndrica

1.0–40.0 mg/L Mo(VI)


1.6–65.0 mg/L MoO₄²⁻




1
4.0 mL Probe
Sample
Echantillon
Muestra




2
Schütteln
Shake
Agiter
Agitar



3
Säubern
Clean
Nettoyer
Limpiar



4
5'00 min



5
Messung
Measurement
Measure
Medición