

Nitrat 50

Nitrate / Nitrato

Test 0-64

REF 985 064



Rundküvettestest / Tube test

Test en cuves rondes / Prueba en cubeta cilíndrica

0.3–22.0 mg/L $\text{NO}_3\text{-N}$


2–100 mg/L NO_3^-

365 nm


Method(e) / Método

0641


0642




500 μL Probe
Sample
Echantillon
Muestra




500 μL R2




Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo



Säubern
Clean
Nettoyer
Limpiar



10'00 min



Messung
Measurement
Mesure
Medición