

# pH 6.5–8.2

# Test 0-72

REF 918 72



Rundküvettestest / Tube test

Test en cuves rondes / Prueba en cubeta cilíndrica


436 / 540 nm

Method(e) / Método


0721

pH 6.5–8.2


Nullwert / Blanc value / Zéro / Cero



10 mL Probe  
Sample  
Echantillon  
Muestra




Säubern  
Clean  
Nettoyer  
Limpiar




Messung  
Measurement  
Mesure  
Medición


Messwert / Sample / Echantillon / Muestra




10 mL Probe  
Sample  
Echantillon  
Muestra




1 x NANOFIX pH




Schütteln  
Shake  
Agiter  
Agitar



Säubern  
Clean  
Nettoyer  
Limpiar



3'00 min



Messung  
Measurement  
Mesure  
Medición