

Sauerstoff 12

Oxygen / Oxygène / Oxígeno

Test 0-82

REF 985 082




Rundküvettentest / Tube test
Test en cuves rondes / Prueba en cubeta cilíndrica

436 nm


Method(e) / Método

0.5–12.0 mg/L O₂


0821




1
Probe
Sample
Echantillon
Muestra




2
Nullwert
Blank value
Zéro
Cero




3
2 x R1




4
2 x R2




5
Schütteln
Shake
Agiter
Agitar




6
2'00 min




7
5 x R3



8
Schütteln
Shake
Agiter
Agitar



9
Säubern
Clean
Nettoyer
Limpiar



10
Messung
Measurement
Mesure
Medición