


NANOCOLOR® Ammonium 3

Amonio



Test 0-03
REF 985003

Methode / Method / Méthode / Método

0031 0.04–2.30 mg/L NH₄-N
0032 0.05–3.00 mg/L NH₄⁺
0033 0.05–3.00 mg/L NH₃

690 nm

Min	20–80 %		Max
0.04	0.49	1.85	2.30 mg/L NH ₄ -N
0.05	0.64	2.41	3.00 mg/L NH ₄ ⁺ / NH ₃

