


# NANOCOLOR® Nitrat 50

## Nitrate / Nitrato



**Test 0-64**  
**REF 985064**

Methode / Method / Méthode / Método  
0641 0.3–22.0 mg/L NO<sub>3</sub>-N  
0642 2–100 mg/L NO<sub>3</sub><sup>-</sup>

345 / 350 / 365 nm

Min	20–80 %		Max
0.3	4.6	17.7	22.0 mg/L NO <sub>3</sub> -N
2	22	80	100 mg/L NO <sub>3</sub> <sup>-</sup>

