


# NANOCOLOR<sup>®</sup> gesamt-Phosphat 15

## total Phosphate / Phosphate total / Fosfato total













**Test 0-80**  
**REF 985080**

Methode / Method / Méthode / Método

0801 0.30–15.00 mg/L P  
0802 1.0–45.0 mg/L PO<sub>4</sub><sup>3-</sup>  
0803 0.7–34.5 mg/L P<sub>2</sub>O<sub>5</sub>


690 nm

Min	20–80 %		Max
0.30	3.24	12.06	15.00 mg/L P
1.0	9.8	36.2	45.0 mg/L PO <sub>4</sub> <sup>3-</sup>



# NANOCOLOR<sup>®</sup> ortho-Phosphat 15

## ortho-Phosphate / ortho Phosphates / orto Fosfato



**Test 0-80**  
**REF 985080**

Methode / Method / Méthode / Método

0805 0.30–15.00 mg/L PO<sub>4</sub>-P

0806 1.0–45.0 mg/L PO<sub>4</sub><sup>3-</sup>

0804 0.7–34.5 mg/L P<sub>2</sub>O<sub>5</sub>

690 nm

Min	20–80 %		Max
0.30	3.24	12.06	15.00 mg/L PO <sub>4</sub> -P
1.0	9.8	36.2	45.0 mg/L PO <sub>4</sub> <sup>3-</sup>

