


# NANOCOLOR® gesamt-Phosphat 5

## total Phosphate/Phosphate total/Fosfato total













**Test 0-81**  
**REF 985081**

Methode / Method / Méthode / Método

0811 0.20–5.00 mg/L P  
0812 0.5–15.0 mg/L PO<sub>4</sub><sup>3-</sup>


690 nm

Min	20–80 %		Max
0.20	1.16	4.04	5.00 mg/L P
0.5	3.4	12.1	15.0 mg/L PO <sub>4</sub> <sup>3-</sup>



# NANOCOLOR<sup>®</sup> ortho-Phosphat 5

## ortho-Phosphate / ortho Phosphates / orto Fosfato



**Test 0-81**  
**REF 985081**

Methode / Method / Méthode / Método

0815 0.20–5.00 mg/L PO<sub>4</sub>-P

0816 0.5–15.0 mg/L PO<sub>4</sub><sup>3-</sup>

690 nm

Min	20–80 %		Max
0.20	1.16	4.04	5.00 mg/L PO <sub>4</sub> -P
0.5	3.4	12.1	15.0 mg/L PO <sub>4</sub> <sup>3-</sup>

