


NANOCOLOR[®] gesamt-Stickstoff TN_b 220

total Nitrogen / Azote total / Nitrógeno total











Test 0-88
REF 985088

Methode / Method / Méthode / Método
0881 5–220 mg/L N


345 / 350 / 365 nm

Min	20–80 %	Max
5	48 177	220 mg/L N



NANOCOLOR[®] gesamt-Stickstoff TN_b 220

total Nitrogen / Azote total / Nitrógeno total



Test 0-88
REF 985088

Methode / Method / Méthode / Método
0881 5–220 mg/L N

345 / 350 / 365 nm

Min	20–80 %	Max
5	48 177	220 mg/L N

