


visocolor® ECO Chrom(VI)

Chromium(VI) / Chrome(VI) / Cromo(VI)




Messbereich / Range / Domaine de mesure / Rango
0.02 - 0.50 mg/l Cr(VI)


REF 931 020
931 220




2 x 5 ml Probe
Sample
Echantillon
Muestra




5 Cr-1




Schütteln
Shake
Agiter
Agitar




5 Cr-2



Schütteln
Shake
Agiter
Agitar



3'00 min



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar