


visocolor® ECO Nitrat

Nitrate / Nitrato




Messbereich / Range / Domaine de mesure / Rango
1 - 120 mg/l NO₃⁻


REF 931 041
931 241



2 x 5 ml Probe
Sample
Echantillon
Muestra




5 NO_3^-1




Schütteln
Shake
Agiter
Agitar




1 NO_3^-2



Kräftig schütteln
Shake well
Bien agiter
Agitar intensamente



5'00 min



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar