

MACHEREY-NAGEL

NANOCOLOR[®] reagent cases

Water Analysis



Mobile photometric analysis

- Modular case solutions for highest flexibility
- Mobile analysis with excellent precision
- The complete equipment in one robust case

MACHEREY-NAGEL

www.mn-net.com



Mobile photometric analysis

Robust reagent cases

The practical **NANOCOLOR**[®] reagent cases have been designed especially for mobile use. These mini labs consisting of stable cases with high-quality foam inlays are available with or without photometers. The high modularity allows equipping the cases with all accessories required for analysis.

In combination with the new **NANOCOLOR**[®] *VARIO Mini*, the reagent cases allow a comfortable and easy analysis on-site. Thus results can be evaluated immediately, even if a preceding sample digestion is necessary.

With the compact photometer **PF-12^{Plus}** the **NANOCOLOR**[®] reagent case becomes an all-rounder. Due to the wide variety of preprogrammed **NANOCOLOR**[®] tube tests and **VISOCOLOR**[®] *ECO* tests on the **PF-12^{Plus}**, the reagent case is perfectly suited for mobile analysis. Equipped with further accessories (heating block, pipettes, etc.) and up to three **NANOCOLOR**[®] tube tests, sample digestions are no longer a barrier on-site.



Mobile photometric analysis

Photometer PF-12^{Plus}

- Mobile compact photometer
- Easy handling for precise results
- Robust and water proof according to IP 68
- Evaluation of *NANOCOLOR*[®] tube tests and *VISOCOLOR*[®] *ECO* tests
- Nephelometric turbidity measurement and NTU-Check
- More than 100 preprogrammed methods

REF 919 250



Photometer PF-3 COD

- Mobile compact photometer
- Usable for mobile COD analysis
- Handy and robust due to compact dimensions
- Easy handling with 4-button-menu navigation
- Robust and water proof according to IP 68

REF 919 342



Heating block *NANOCOLOR*[®] *VARIO C2*

- User-friendly touch screen with intuitive icons
- Extremely short heating time
- Simultaneous digestion of up to 12 samples
- High temperature stability
- Graphically displayed heating curves
- Internal quality control with the *NANOCOLOR*[®] T-Set or *NANOCOLOR*[®] USB T-Set

REF 919 350



Heating block *NANOCOLOR*[®] *VARIO Mini*

- Flexible power supply for mobile sample preparation
- User friendly handling with 4-button-menu navigation and intuitive icons
- Simultaneous digestion of up to 6 samples
- Reliable temperature stability
- High quality heating block with perfect price-performance ratio
- Internal quality control with the *NANOCOLOR*[®] USB T-Set

REF 919 380



Mobile photometric analysis

NANOCOLOR® VARIO Mini

The NANOCOLOR® VARIO Mini is particularly suitable for mobile use. Due to its 4-button-menu navigation, working is easy and comfortable. In combination with our photometers and the appropriate NANOCOLOR® tube test, a fast and reliable information on the water quality can be obtained on-site. To ensure high flexibility with our system, the NANOCOLOR® VARIO Mini is available separately in a cardboard box. Together with the NANOCOLOR® reagent case, the NANOCOLOR® VARIO Mini is adaptable even to specific customer requirements.

Your benefit

- Handy and compact heating block for mobile analysis
- Digestion of up to 6 samples simultaneously
- NANOCOLOR® USB T-Set for temperature control and calibration of the heating block

Ordering information

Description	REF
NANOCOLOR® reagent case with photometer PF-12 ^{plus} with additional space for 1 heating block NANOCOLOR® VARIO C2 or NANOCOLOR® VARIO Mini, 2 pipettes and 3 NANOCOLOR® tube tests	919 214
NANOCOLOR® reagent case PF-12 ^{plus} , empty with additional space for 1 photometer PF-12 ^{plus} , 1 heating block NANOCOLOR® VARIO C2 or NANOCOLOR® VARIO Mini, 2 pipettes and 3 NANOCOLOR® tube tests	919 215
NANOCOLOR® reagent case with photometer PF-3 COD with additional space for 1 heating block NANOCOLOR® VARIO C2 or NANOCOLOR® VARIO Mini, 2 pipettes and 3 NANOCOLOR® tube tests	919 212
NANOCOLOR® reagent case PF-3 COD, empty with additional space for 1 photometer PF-3 COD, 1 heating block NANOCOLOR® VARIO C2 or NANOCOLOR® VARIO Mini, 2 pipettes and 3 NANOCOLOR® tube tests	919 213



NANOCOLOR® VARIO Mini



For the power supply the NANOCOLOR® VARIO Mini can be easily connected to the voltage of the car's electrical system.



Example: Fully assembled NANOCOLOR® reagent case with compact photometer PF-12^{plus}



Example: Fully assembled NANOCOLOR® reagent case with compact photometer PF-3 COD

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

