

Filter Papers – Filtrierpapiere

Grade / Sorte	MN 214ff
Weight / Gewicht [g/m²]	140
Thickness / Dicke [mm]	0,28
Mirgration distance / Saughöhe [mm/10 min]	55 - 60
Surface / Oberfläche	smooth / glatt
Applications and properties /	For determination of fats, fat-free, low nitrogen content, low ash content, for paper chromatography /
Anwendungen und Eigenschaften	Zur Fettbestimmung, fettfrei, stickstoffarm, aschearm für Papierchromatographie

These data are typical values and must not be regarded as specification.
Die angegebenen Daten sind typische Werte und stellen keine Spezifikation dar.